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

Acculis


Microwave Tissue Ablation System

A single device for fast,
predictable, large ablations.


 **angiodynamics**

Considerations in selecting a Microwave Ablation System


Simple set up, planning, device placement and coagulation parameters

Planning


Placement

CT


Post-ablation

Ex-vivo tissue, bovine

A single device for fast, predictable, large ablations.

Fast


The Acculis* MWA system is rated the fastest in achieving the largest coagulation diameter (short-axis) at 5 and 6 minutes.¹

Line graphs show the (a) volume and (b) short-axis diameter of the coagulation zone in relation to ablation duration.

• Acculis (@ 180 watts)	• $b = 17.4 (SE 1.4) * \ln(t) + 13.1 (SE 2.2)$
• HS Amica	• $b = 11.7 (SE 0.6) * \ln(t) + 18.3 (SE 1.4)$
• Covidien, 1 Antenna	• $b = 8.0 (SE 0.6) * \ln(t) + 8.5 (SE 1.0)$
• Covidien, 3 Antennas	
• Medwaves	• $b = 8.9 (SE 0.8) * \ln(t) + 11.5 (SE 1.8)$

Ex-vivo bovine study—liver

Predictable

The Acculis MWA system provides the most spherical and largest diameter ablations at 5 minutes.^{1,2}

Graphs show the volume, b/a ratio and short axis diameter of the coagulation zone after 5 minutes, 10 minutes, and the maximum recommended ablation duration:

System A — Acculis, 6 minutes
 System B — HS Amica, 10 minutes
 System C — Covidien, 10 minutes
 System D — Medwaves, 15 minutes

1 Hoffman et al., Radiology, Published online, Feb 2013
 2 IFU pMTA Issue 4
 3 IFU MTA Issue 4

Acculis MTA System

- Thermal coagulation of Soft Tissue using Microwave Energy
- Probe placed under CT or Ultrasound Guidance
- System consists of the Sulis VpMTA Generator, Acculis Local Control Station (LCS), Accu2i pMTA Applicators, and optional MTA Temperature Probes.

Accu2i pMTA Applicator

Percutaneous applicator (available in 14 cm, 19 cm and 29 cm lengths)


- Fast, large, near-spherical ablations¹
- 4.5 cm x 5.5 cm diameter ablations in 6 min²

Power	2 minutes	4 minutes	6 minutes
Liver			
60 watts	2.3 x 2.7	3.3 x 3.7	3.5 x 4.0
100 watts	2.7 x 3.2	3.6 x 4.1	4.0 x 4.7
140 watts	2.9 x 3.5	4.1 x 5.2	4.5 x 5.5

Diameter x Length (both in cm)

See IFU for Muscle and Kidney Protocols with the applicators.

NOTE: Coagulation results are from coagulations performed in ex-vivo animal tissue models

Accu20s MTA Surface Applicator

Open use (available in 30 cm length)

- 1.9 cm x 1.2 cm diameter ablations in 60 sec²
- Scalable burn size³
- Unique design—only microwave surface applicator on the market

Power	20 seconds	40 seconds	60 seconds
80 watts	0.8 x 0.6	1.4 x 0.9	1.6 x 1.0
90 watts	1.2 x 0.6	1.5 x 0.9	1.7 x 1.1
100 watts	1.5 x 0.8	1.7 x 0.9	1.9 x 1.2

Diameter x Depth (both in cm)

NOTE: Coagulation results are from coagulations performed in ex-vivo animal tissue models

Accu5i MTA Applicator

Open or Laparoscopic use (available in 17.4 cm and 33.4 cm lengths)


- 5.0 cm x 6.8 cm diameter ablations in 8 min³
- 5.6 mm shaft³
- Scalable near spherical coagulation zones³

Power	4 minutes	6 minutes	8 minutes
Liver			
30 watts	2.4 x 2.8	2.9 x 3.3	3.6 x 3.9
70 watts	3.4 x 4.8	4.0 x 5.5	4.7 x 6.3
100 watts	4.3 x 6.0	4.6 x 5.9	5.0 x 6.8

Diameter x Length (both in cm)


See IFU for Muscle and Kidney Protocols with the applicators.

NOTE: Coagulation results are from coagulations performed in ex-vivo animal tissue models


Acculis LCS (Local Control Station)

- Coolant temperature monitoring system
- Most efficient pump on the market
- Integrated control interface for disposable applicators
- Single-action multi-connector cartridge

Sulis V^{pMTA} Generator

- Touch screen user interface
- Most powerful generator on the market
- Simple to operate¹
- 140 watts maximum

Acculis MTA temperature probes

- Single-use, disposable instrument used to monitor temperature levels up to 70 °C in real-time


Acculis

Microwave Tissue Ablation System

A single device for fast,
predictable, large ablations.

4.5 cm x 5.5 cm ablation[†] in a single
6 minute application

[†]The liver coagulation results are from coagulations
performed in ex-vivo animal tissue models.

Product Name	Diameter	Length	Part #
pMTA Accu2i Standard Applicator (cooled)	1.8 mm	14 cm	900-600
pMTA Accu2i Intermediate Applicator (cooled)	1.8 mm	19 cm	900-602
pMTA Accu2i Long Applicator (cooled)	1.8 mm	29 cm	900-601
MTA Accu5i, Applicator (not cooled)	0.56 cm	17.4 cm	900-312
MTA Accu5i, Applicator (not cooled)	0.56 cm	33.4 cm	900-313
MTA Accu20s, Surface Applicator (not cooled)	24 mm	30 cm	900-500
MTA Temperature Probes (black and grey)	18G	20 cm	900-314
MTA Temperature Probes (black and grey)	18G	33 cm	900-315
Sulis V ^{pMTA} Generator	—	—	806-010
LCS Unit	—	—	806-011
Cart and Monitor Arm	—	—	700-006
Cord Set Canada 10A/125v	—	—	800-017

IMPORTANT RISK INFORMATION

DFU (510-243)

The Accu2i pMTA applicator is intended for the coagulation of soft tissue during open, laparoscopic or percutaneous procedures.

DFU (62833-09-002-20s) & (62810-09-002-Accu5i)

The Acculis Accu20s MTA Surface Applicator and Accu5i MTA Surgical Applicator are indicated for the intraoperative coagulation of soft tissue.

CONTRAINDICATIONS: The Accu2i pMTA applicators, Acculis Accu20s MTA Surface Applicators and Accu5i MTA Surgical Applicators are contraindicated in patients with heart pacemakers and other electrical device implants.

For exclusive international use, outside the US.

International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands
tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

www.angiodynamics.com

*AngioDynamics, the AngioDynamics logo, Acculis and Sulis are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary.

Legal Manufacturer:
AngioDynamics, Inc.
Parklands Business Park
Denmead, Hampshire UK

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We look forward to your favourable consideration and assure you of our best service at all time.

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
CE
0123



Microwave Ablation A New Portable Generation

About Canyon Medical


Canyon Medical is a Chinese company based in Nanjing since 1986, where we manufacture a range of microwave medical equipment of our own design. We are specialized in Microwave Ablation for about 20 years with more than 200,000 accumulated cases done successfully with our system in China. By innovating with clinical practices hands in hands, we are pioneering to solve clinical challenges in this field. In view of our continuous and prominent contribution, our company was granted with many great honors, the most prominent award is National Award for Scientific and Technological Invention in 2014.

In 2015, Canyon Medical joined the family of Micro-Tech. In virtue of Micro-Tech's international platform and management philosophy, Canyon Medical is advancing towards a greater international stage at a higher level.


As a leading brand in China with endless efforts in product design and R&D, we are concentrating on the development of high value-added microwave ablation products that meets customers' needs.


Chinese Medical Science and Technology Prize, Second Prize, 2012

Chinese Medical Science and Technology Prize, Second Prize, 2018


National Technology Invention Award, Second Price, 2014

Canyon as the exclusive sponsor of Chinese Microwave Ablation Master Summit since 2015

MICROWAVE ABLATION GENERATOR


- A new portable solid state system to ensure stable microwave output.
- User friendly Interactive and sensitive touchscreen with adjustable power, timing and temperature setting.
- Continuous and pulsed modes are available for different clinical options.
- Incorporated multi-safety monitoring design.

Ordering and Specifications


Microwave Ablation Generator					
REF.	Frequency	Display	Energy Output	Time Setting	Working Power
KY-2000A	2450 MHz	7 in. Touchscreen	5-100 W	1-30 Minutes	100-240 V 50-60 Hz

MICROWAVE ABLATION ANTENNAS


- High performance material to ensure better connection and strength.
- Unique shaft coating material to avoid tissue adhesion.
- Sharp tip needle and thin diameter, less invasive.
- Internally water cooled antenna with intelligent temperature control.
- Signal enhancement design for better visibility under CT or ultrasound guidance.
- High efficiency detachable cable, unique enhanced connection with the antenna.

THYTIP™ Thyroid series

REF	Diameter	Active Tip	Shaft Length
KY-2450B-3	17 G	3 mm	100 mm
KY-2450B-4	18 G	3 mm	70 mm
KY-2450B-5	16 G	3 mm	100 mm
KY-2450B-13	16 G	3 mm	70 mm

SURETIP® Standard series

REF	Diameter	Active Tip	Shaft Length
KY-2450B-T1	16 G	11 mm	100 mm
KY-2450B-T2	16 G	11 mm	150 mm
KY-2450B-T3	16 G	11 mm	200 mm
KY-2450B-T4	15 G	11 mm	100 mm
KY-2450B-T6	15 G	11 mm	150 mm
KY-2450B-T7	15 G	11 mm	200 mm
KY-2450B-T8	15 G	11 mm	180 mm
KY-2450B-T9	14 G	11 mm	180 mm
KY-2450B-T10	15 G	11 mm	250 mm
KY-2450B-T5	18 G	11 mm	150 mm
KY-2450B-T11	18 G	11 mm	100 mm
KY-2450B-T13	18 G	11 mm	200 mm
KY-2450A-7	16 G	5 mm	100 mm
KY-2450A-11	16 G	5 mm	130 mm
KY-2450B-12	15 G	5 mm	180 mm


XPHRETIP® Spherical series


REF	Diameter	Active Tip	Shaft Length
KY-2450B-QT	14 G	25 mm	200 mm
KY-2450B-QT1	14 G	25 mm	100 mm
KY-2450B-QT2	14 G	25 mm	150 mm
KY-2450B-QT3	14 G	25 mm	250 mm
KY-2450B-QT4	15 G	25 mm	100 mm
KY-2450B-QT5	15 G	25 mm	150 mm
KY-2450B-QT6	15 G	25 mm	200 mm
KY-2450B-QT7	15 G	25 mm	250 mm
KY-2450B-QT8	15 G	25 mm	300 mm

over 20 years with 200,000+ accumulated cases

Precise Stable Safer


40w


5min 2.5cm*2.9cm*0.1cm
10min 3.0cm*3.9cm*0.2cm
15min 3.3cm*4.3cm*0.3cm

50w


5min 2.7cm*3.2cm*0.2cm
10min 3.1cm*4.4cm*0.2cm
15min 3.7cm*4.5cm*0.3cm

60w


5min 2.7cm*3.5cm*0.2cm
10min 3.5cm*4.0cm*0.2cm
15min 3.8cm*4.6cm*0.3cm

70w


5min 2.8cm*4.3cm*0.3cm
10min 3.5cm*4.8cm*0.3cm
15min 4.0cm*5.1cm*0.4cm

80w


5min 3.0cm*4.9cm*0.3cm
10min 3.6cm*4.9cm*0.4cm
15min 4.3cm*5.2cm*0.5cm

100w

5min 3.2cm*4.9cm*0.4cm
10min 3.9cm*5.0cm*0.4cm
15min 4.4cm*5.6cm*0.6cm


Ablation by two antennas working simultaneously

Note: Coagulative performance by 15G MW antenna carried on ex-vivo porcine liver.
The actual ablation dimensions in in-vivo tissues may differ from such static laboratory results.

BENEFITS OF MICROWAVE ABLATION

Microwave ablation can generate larger ablation zone in shorter time with less susceptibility to heat sink effect. Moreover, MWA radiates all biological issues, such as bone, lung, charred or desiccated tissues. Besides, multiple antennas can be placed to treat simultaneously for large tumors or multifocal lesion.

Furthermore, with Canyon New Generation Microwave Ablation, the achieved ablation zone can be uniform, homogeneous and predicable.

MWA in a 45-year-old woman with HCC on right lobe of the liver.

Notes: (a) Preablation conventional ultrasound (US) scan shows a hypoechoic lesion with blood supply (arrows). (b) Contrast-enhanced US before MWA shows tumor enhancement in arterial phase with the size of 5.5×4.7 cm (arrows). (c) Conventional US shows two microwave antennas (mark) are placed in the tumor. (d) Contrast-enhanced US shows no enhancement of the ablation zone at 12 month after treatment (arrows). (e) MRI scan obtained 21 month after ablation shows hypoattenuating ablation zone (arrow) without enhancement.^①

① P. Liang et al. (eds.), Microwave Ablation Treatment of Solid Tumors, DOI 10.1007/978-94-017-9315-5_4,
© Springer Science+Business Media Dordrecht 2015

CANYON MEDICAL INC.

Add: Building 3, Phase 2 Accelerator, No. 11 Yaogu Avenue,
Jiangbei New Area, 210032, Nanjing, Jiangsu, China
Email: sales@canyonmedical.com.cn
Website: www.canyonmedical.com.cn

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